

PRODUCT	DESCRIPTION	VISCOSITY (cP)*	THIX INDEX	GEL TIME (min)	Promotor / Accelerator (% m/m)	Organic Peroxide Catalyst (% m/m)	HDT (°C)
OrthoPhthalic / General Purpose Laminating Resins							
UPPOL-3404TA-H30	Ortho, Fast Wet-out, medium reactive Thixotropic Laminating	500-600	1.3-1.6	28-32	PA	1.0	62
UPPOL-3877TA H30	Terephthalic fast Wet-out, Thixotropic Laminating	400-500	1.3-1.5	30	PA	1.0	82
UPPOL-3401TA-H30	Ortho, Medium Reactive Thixotropic Laminating	850-950	1.1-1.3	30	PA	1.0	65
UPPOL-3562-230	Ortho, Medium Reactive, unaccelerated, Non Thix, Laminating	250-350	-	14-22	1% (1% Co Oct)	1.0	86

*Refer to Technical Data sheet for methods

PA = Pre-Accelerated

HDT: Method ISO 0075-B